



BORT MEDICAL

has developed a

HIGH QUALITY BANDAGE

with interchangeable

Hinges and Splints

Product Description

- * Patella fixation bandage
- * Anatomically form-knitted
- * Strong Dual Tension quality
- * Pressure-reduced areas at edges and at hollow of the knee
- * With silicone pad
- * One Removable Hinge each side
- * PLUS One splint each side.
- * Hinges Provide Extension limit at 0°.
- * Hinges and Splints are Interchangeable as needed.
- * Excellent fit thanks to innovative plus size system.

Indications

Patella instability, articular effusion, swellings, arthrosis and arthritis, post-operative irritations, ligament instability, gonarthrosis, Lateral ligament instability, chondropathia patellae.

Sizes

d-measure	f-measure
26 - 29 cm	32 - 36 cm- small
29 - 33 cm	36 - 40 cm- medium
29 - 33 cm	40 - 44 cm- medium plus
33 - 37 cm	40 - 44 cm- large
33 - 37 cm	44 - 48 cm- large plus
37 - 41 cm	44 - 48 cm- x-large
37 - 41 cm	48 - 52 cm- x-large plus
41 - 45 cm	52 - 56 cm- xx-large
41 - 45 cm	56 - 60 cm- xx-large plus
45 - 49 cm	60 - 64 cm- xxx-large
45 - 49 cm	64 - 68 cm- xxx-large plus

Color: Silver

1 cm = 0,3937 inch

What's NEW at BORT MEDICAL?



StabiloGen® Eco with Articulated Joint/ Splint Combo Article No.: 114 520J

L1820

*NEWEST ADDITION: THIGH STRAP for increased Stability

BSOS

Bort-Swiss Orthopedic Supply, LLC
2221 Applebrook Dr.
Commerce MI 48382
248-242-6907-Office
248-242-6908-Fax

beth.swiss@bort-swissortho.com
www.bort-swissortho.com

Locate us on these websites:
www.bort-swissortho.com
www.facebook.com
www.twitter.com

The complete Bort Medical line can be viewed on the Bort website.
www.bort.com