



NOVELTIES 2006



... FOR SHOULDER, ARM, BACK,
KNEE AND ANKLE



NOVELTIES 2006

IN COOPERATION WITH DOCTORS AND ORTHOPAEDIC TECHNICIANS ALREADY EXISTING PRODUCTS COULD BE IMPROVED AND NEW PRODUCTS HAVE BEEN DEVELOPED:



NEW

Art. no. 112 510

BORT MobiDig® Extension

- > Finger splint for extension of PIP joints
- > Newly designed sizes system for increased individual care
- > Linear spring force
- > Spring force can be adjusted individually
- > High wearing comfort by elaborated, wear-resistant, foamed padding
- > Color: white



NEW

Art. no. 121 100

BORT OmoFX

- > Reusable shoulder-arm abduction orthosis
- > Pleasant smooth velour cover
- > Physiological strap system, thus even increased wearing comfort
- > Immobilization of shoulder joint and elbow
- > Abduction position: approx. 30 °
- > Color: blue



NEW

Art. no. 122 500

BORT EpiContur® Epicondylitis-Splint with 2 Pelottes

- > Special fit for friction or compression of tendon insertion points
- > Two freely positionable pelottes
- > High wearing comfort due to skin-friendly materials
- > Individual pressure-setting by means of infinitely variable adjustable reversing strap
- > Colors: blue, beige



NEW

Art. no. 104 710

BORT Stabilo® Back Bandage with Pelotte

- > Anatomic form-knitted bandage
- > With lumbosacral friction pelotte
- > 2 anatomically pre-cast splints in dorsal area
- > The semi-elastic abdominal position provides correction of the hollow of back of LWS (lower lumbar spine pain)
- > Easy putting on by means of hand loops
- > Colors: beige, silver-grey



NEW

Art. no. 101 600

BORT OA® Spirit GII Unloader™

- > Knee-Unloader orthosis in case of isolated lateral/medial gonarthrosis
- > Varus-/Valgus-correction to be freely adjusted
- > Thermoplastic, easy adjustable thigh and lower leg shells
- > Adjustable, polycentric gear segment hinge
- > Low own weight
- > Individual dosage of axis correction by DNS (Dynamic-Force-Strap)
- > Color: black



NEW

Art. no. 054 900

BORT AchilloStabil® Ankle Bandage

- > Three-dimensional form-knitted dual tension bandage
- > Visco-elastic profile inserts alongside of the achilles tendon
- > Flat-knitting quality with specially padded seams
- > Pressure relief of achilles tendon and intermittent massage for better blood circulation
- > With heel cushion for non-effected leg to equalize length differences of legs
- > Colors: beige, silver-grey

QUALITY
MANAGEMENT
certified acc. to
medical products guidelines
93/AZ/EWG



QUALITY
MANAGEMENT
we are certified
regular re-auditing
certified acc. to ISO 9001:2000



BORT GmbH
P.O. Box 1330
D-71367 Weinstadt
medical@bort.de
www.bort.de

Phone **49 7151 99200-0
Fax **49 7151 99200-50