


In order to manipulate gripping more easily

Mobility is not always given! Many people suffer from restricted mobility. BORT has developed a special reacher to facilitate daily activities, so that no insurmountable problem arises for them. As independence and safety should be a matter-of-course.



Bort Reachers are awarded with the „Red Point“ of the Design Center in Essen for outstanding design.

Ideal for anywhere ►
and every situation:
The 75 cm long foldable
version can be collapsed
to 40 cm, is space-saving
and very handy.



▲ The BORT Reacher re-allows patients with restricted mobility and senior people to easily manage daily activities again. Gripping of small objects or picking up dropped items are no longer a problem with the BORT Reacher!



Safe gripping:



The new type of reachers for patients with restricted mobility

Even of smallest items. No slipping by extra large contact surface. Solid grip thanks to locking device of the lever handle. Simple, controlled compression.

Two different adjustment possibilities:



Short leverage: normal effort for unrestricted manual gripping function.
Extended leverage: reduced effort for restricted manual gripping function.

The useful locking device:



Grips items safely even when the handle is released and is integrated in all versions.

Always along:



The 75 cm long foldable version can be collapsed to 40 cm.

For every situation:



Order-Information:	Art. no.	Length
Foldable	118 100	75/40 cm
Standard	118 150	50 cm
	118 170	75 cm
	118 190	90 cm

YOUR SPECIALIZED DEALER:

249-CreitmA/GB-04/03.wm



BORT GmbH
 P. O. Box 1330
 D-71367 Weinstadt-Benzach
 www.bort.de



Simply grip it

